

Product Specification Sheet

Anti-Protein Phosphatase 1 γ 1
Product Number: PP15
Aliquot: 100 μ g
Lot Number pp15.10552A
Storage: -70°C

DESCRIPTION:	This is a rabbit polyclonal antibody is specific for serine/threonine protein phosphatase 1 γ 1 (PP 1 γ 1) corresponding to the C-terminus region. This antibody does not cross-react with other protein phosphatase isoforms.
MOLECULAR WEIGHT:	37 kDa
CONCENTRATION:	1 mg/mL determined by BCA method using BSA as a standard.
STORAGE BUFFER:	Phosphate Buffered Saline; 0.08% (w/v) sodium azide; Sterile filtered solution.
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	NH ₂ -Thr-Pro-Pro-Arg-Gly-Met-Ile-Thr-Lys-Gln-Ala-Lys-Lys-COOH Conjugated to KLH
CLASS:	IgG
PURIFICATION:	Ammonium Sulfate Fractionation
APPLICATIONS:	This antibody is suitable for use in Western blotting at 5-10 μ g/mL. Reported utility in immunoprecipitation.
SPECIES SPECIFICITY:	Human, mouse, rat, bovine.