

Product Specification Sheet

Anti-Serine/Threonine Protein
Phosphatase 2 A/A
Product Number: PP19
Aliquot: 100 µg
Lot Number pp19.8412
Storage: -70°C

DESCRIPTION:	This is a rabbit polyclonal antibody specific for serine/threonine protein phosphatase 2 A/A corresponding to the N-terminal peptide of PP2A/A regulatory subunit. This antibody does not cross-react with other protein phosphatase isoforms.
MOLECULAR WEIGHT:	65 kDa
CONCENTRATION:	1 mg/mL determined by BCA method using BSA as a standard.
STORAGE BUFFER:	Phosphate Buffered Saline; 0.08% (w/v) sodium azide; Sterile filtered solution.
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	AA sequence 7-19.
CLASS:	IgG
PURIFICATION:	IgG Purified.
APPLICATIONS:	This antibody is suitable for use in Western blotting at 5-10 µg/mL and immunocytochemistry at 10-15 µg/mL.
SPECIES SPECIFICITY:	Human, mouse, rat, bovine.