

Product Specification Sheet

Anti-Protein Phosphatase 2
A/B γ 2
Product Number: PP 21
Aliquot: 100 μ g
Lot Number pp21.10435
Storage: -70°C

| | |
|-----------------------------|---|
| DESCRIPTION: | This is a rabbit polyclonal antibody specific for serine/threonine protein phosphatase 2 A/B γ 2 (PP2A/B γ 2) corresponding to amino acids 53-66 of PP2A/B. This antibody does not cross-react with other protein phosphatase isoforms. |
| MOLECULAR WEIGHT: | ~50 kDa |
| CONCENTRATION: | 1 mg/mL determined by BCA method using BSA as a standard. |
| STORAGE BUFFER: | Phosphate Buffered Saline; 0.08% (w/v) sodium azide; Sterile filtered solution. |
| STORAGE: | -70°C. AVOID REPEATED FREEZE THAWS. |
| IMMUNOGEN: | Amino acids 53-66 of PP2A/B |
| ISOTYPE: | IgG |
| PURIFICATION: | Ammonium Sulfate Fractionation |
| APPLICATIONS: | This antibody is suitable for use in Western blotting at 5-10 μ g/mL. |
| SPECIES SPECIFICITY: | Human, mouse, rat, bovine. |