

**Product Specification Sheet**

Description: Anti Rabbit P4502E1  
PAB  
Product Number: PR32  
Aliquot: 100uL  
Lot Number: pr32.031996SK  
Storage: -20C

**DESCRIPTION:** Polyclonal goat IgG antibody against rabbit P4502E1

**MOLECULAR WEIGHT:** 57kDa

**RECOMMENDED DILUTION:** 1:1000 for WB

**STORAGE:** -20C

**IMMUNOGEN:** Rabbit P4502E1

**PURIFICATION:** IgG purification

**CLASS:** IgG

**APPLICATIONS:** Western blotting, ELISA

**SPECIES SPECIFICITY:** Rabbit, Rat