

Product Specification Sheet

Anti-Human PAI-1 MoAB
Product Number: PT10
Aliquot: 1.0 mg
Lot Number pt10.521
Storage: -70°C

DESCRIPTION:	This is a non-inhibitory mouse monoclonal antibody against human plasminogen activator inhibitor-1 (PAI-1).
MOLECULAR WEIGHT:	43 kDa non-glycosylated
CLONE:	671-3E5-29
CONCENTRATION:	1.3 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Recombinant full length wild-type human PAI-1.
CLASS:	IgG
PURIFICATION:	Purified by Immobilized Protein G.
APPLICATIONS:	This antibody is suitable for use in Western blotting and ELISA.
SPECIES SPECIFICITY:	Human, Rabbit and Pig