

Anti-Mouse tPA IgG Fraction
Product Number: PT31
Aliquot: 1.0 mg
Lot Number pt31.404A
Storage: -70°C

DESCRIPTION: This polyclonal antibody is produced in rabbit and is presented as an IgG preparation.

MOLECULAR WEIGHT: 70 kDa

CONCENTRATION: 3.9 mg/mL determined spectrophotometrically.

STORAGE BUFFER: 0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6

STORAGE: -70°C. **AVOID REPEATED FREEZE THAWS.**


IMMUNOGEN: Recombinant Mouse tPA

CLASS: IgG

PURIFICATION: Purified by Protein A Sepharose.

APPLICATIONS: Western: 1-5 µg/mL

SPECIES SPECIFICITY: Mouse.



Non-reduced:
-Lane 1: PT 31 (0.5 µg)
-Lane 2: PT 31 (0.5 µg) + PI 05 (3 µg)
Reduced:
-Lane 3: PT 31 (0.5 µg)
-Lane 4: PT 31 (0.5 µg) + PI 05 (3 µg)
-Lane 5: MW ladder