

**Product Specification Sheet**

Anti-Mouse tPA Biotin Labeled  
IgG  
Product Number: PT32  
Aliquot: 1.0 mg  
Lot Number pt32.1011  
Storage: -70°C

<b>DESCRIPTION:</b>	This polyclonal antibody is produced in rabbit and is biotin labeled.
<b>MOLECULAR WEIGHT:</b>	70 kDa
<b>CONCENTRATION:</b>	3.4 mg/mL determined spectrophotometrically.
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
<b>STORAGE:</b>	-70°C. <b>AVOID REPEATED FREEZE THAWS.</b>
<b>IMMUNOGEN:</b>	Recombinant Mouse tPA
<b>CLASS:</b>	IgG
<b>PURIFICATION:</b>	Purified by Protein A Sepharose
<b>APPLICATIONS:</b>	Western: 5 µg/mL
<b>SPECIES SPECIFICITY:</b>	Mouse