



## Product Specification Sheet

Anti-Mouse tPA Affinity Purified  
High Titer  
Product Number: PT34  
Aliquot: 0.1 mg  
Lot Number pt34.1011b  
Storage: -70°C

<b>DESCRIPTION:</b>	This polyclonal antibody is produced in rabbit and has been affinity purified to produce this high titer aliquot.
<b>MOLECULAR WEIGHT:</b>	70 kDa
<b>CONCENTRATION:</b>	1.7 mg/mL determined spectrophotometrically.
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
<b>STORAGE:</b>	-70°C. <b>AVOID REPEATED FREEZE THAWS.</b>
<b>IMMUNOGEN:</b>	Recombinant Mouse tPA
<b>CLASS:</b>	IgG
<b>PURIFICATION:</b>	Affinity purified by immobilized tPA
<b>APPLICATIONS:</b>	Western: 1:10,000 ELISA: 1:100,000
<b>SPECIES SPECIFICITY:</b>	Mouse.