

Product Specification Sheet

Description: Anti-Human VLDL Receptor Monoclonal Antibody
Product Number: RA32
Aliquot: 1.0 mg
Lot Number: ra32.323A
Storage: -70°C

DESCRIPTION:	Anti-Human VLDL Receptor Monoclonal Antibody. Binds to an epitope on the C-terminus of VLDL Receptor (AA 191-355). Blocks apoE4 binding.
MOLECULAR WEIGHT:	85 kDa
CONCENTRATION:	1.22 mg/mL
STORAGE BUFFER:	0.05M Sodium Phosphate, 0.1M NaCl, 1mM EDTA, pH 6.6
STORAGE:	-70°C
IMMUNOGEN:	Human VLDL Receptor
PURIFICATION:	Immobilized protein A
CLASS:	IgG1 kappa
APPLICATIONS:	ELISA, Western Blotting (non-reducing), Blocks apoE4 binding.
SPECIES SPECIFICITY:	Human. Cross reactivity to other species has not been tested.