

Product Specification Sheet

Description: Recombinant
Porcine tPA, Active
Product Number: UP59
Aliquot: 0.1 mg
Lot Number: up59.413
Storage: -70°C

DESCRIPTION:	This active porcine tPA was recombinantly produced in a stable insect cell line.
MOLECULAR WEIGHT:	70 kDA
CONCENTRATION:	1.15 mg/mL
STORAGE BUFFER:	0.4 M HEPES; 0.1 M NaCl; pH 7.4
STORAGE:	Store at -70°C
PURIFICATION:	Purified by immune affinity chromatography and sterile filtered
SPECIFIC ACTIVITY:	900,000 IU/MG
PURITY:	> 95% by SDS-PAGE analysis
APPLICATIONS:	ELISA, WB