



Product Specification Sheet

Description: Sheep Anti-Mouse and Rat tPA, IgG Fraction, FITC Labeled

Product Number: UP64

Aliquot: 1.0 mg

Lot Number: up64.619

Storage: -70°C

DESCRIPTION:	Sheep Anti-Mouse and Rat tPA, IgG Fraction, FITC Labeled
MOLECULAR WEIGHT:	Target MW is 70 kDa
CONCENTRATION:	8.34 mg/ml
LABEL:	2.2 moles FITC per mole antibody
STORAGE BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
STORAGE:	Store at -70C and protect from light
IMMUNOGEN:	Recombinant full length, glycosylated, wt mouse protein
PURIFICATION:	IgG fraction purified by immobilized Protein G
CLASS:	IgG
HOST:	Sheep
APPLICATIONS:	ELISA
SPECIES SPECIFICITY:	Mouse, Rat