



Product Specification Sheet

RAW Ram Seminal Vesicles
(RSV)
Product Number: V01
Packaging: 250 g
Lot Number: v01.090930
Storage: -70°C

DESCRIPTION:

Ram Seminal Vesicles are a strong source for eicosanoid metabolizing enzymes such as cyclooxygenase, prostaglandins, and numerous other enzymes. These are presented as raw tissue, excised from male ovine bucks immediately after slaughter and stored at -70 °C.

STORAGE:

-70°C **AVOID MULTIPLE FREEZE-THAW CYCLES.**

APPLICATIONS:

This raw prep of RSV's has been used as a resource for eicosanoid metabolizing enzymes.