

**Product Specification Sheet**

Description: Sheep Anti-mouse  
uPAR IgG Fraction, Biotin  
Labeled  
Product Number: VA66  
Aliquot: 1.0 mg  
Lot Number: va66.512  
Storage: -70°C

<b>DESCRIPTION:</b>	Sheep anti-mouse soluble uPAR polyclonal antibody. Biotin labeled.
<b>TARGET:</b>	Recognizes soluble urokinase receptor (sUPAR)
<b>TARGET MOLECULAR WEIGHT:</b>	33 kDa
<b>CONCENTRATION:</b>	10.0 mg/mL
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
<b>STORAGE:</b>	Store at -70°C
<b>IMMUNOGEN:</b>	Recombinant full length mouse uPAR histidine tagged at the N-terminal end
<b>PURIFICATION:</b>	Purified by immobilized protein G
<b>CLASS:</b>	IgG
<b>APPLICATIONS:</b>	ELISA, WB
<b>SPECIES SPECIFICITY:</b>	Mouse