

Product Specification Sheet

Description: Rat Factor X
Product Number: CF42
Aliquot: 0.1 mg
Lot Number: Typical Lot
Storage: -20°C

DESCRIPTION:	Rat Factor X
MOLECULAR WEIGHT:	63 kDa
CONCENTRATION:	5.0 mg/mL
STORAGE BUFFER:	50% glycerol
STORAGE:	Store at -20°C
PURIFICATION:	Prepared from fresh frozen plasma and purified by immunoaffinity chromatography.
PURITY:	> 95% by SDS PAGE analysis