

Product Specification Sheet

Description: Sheep Anti-Human
Factor XII IgG Fraction, Biotin
Labeled

Product Number: CF49

Aliquot: 1.0 mg

Lot Number: cf49.1112

Storage: -70°C

DESCRIPTION:	Sheep Polyclonal Anti-Human Factor XII IgG Fraction, Biotin Labeled
TARGET MOLECULAR WEIGHT:	80 kDa
CONCENTRATION:	5.15 mg/ml
STORAGE BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
STORAGE:	Store at -70°C
IMMUNOGEN:	Full length native human protein
PURIFICATION:	IgG Fraction
APPLICATIONS:	ELISA, WB
SPECIES SPECIFICITY:	Human