

Product Specification Sheet

Description: Anti-human TGF
beta 1 monoclonal antibody
Product Number: CK50
Aliquot: 1.0 mg
Lot Number: Typical Lot
Storage: -20°C

DESCRIPTION: Monoclonal antibody against Human Transforming Growth Factor Beta 1 (TGF-beta 1)

MOLECULAR WEIGHT: 25 kDa

STORAGE BUFFER: 0.5 x PBS, pH 7.2

STORAGE: Once reconstituted, aliquot and store at -80°C.

IMMUNOGEN: Recombinant Human TGF-beta 1.

PURIFICATION: IgG1 antibody was purified by Protein G affinity chromatography.

CLASS: IgG1

PURITY: >98%

APPLICATIONS: ELISA 0.2 ng/well.

SPECIES SPECIFICITY: Human.