

Product Specification Sheet

Description: Anti-human BMP-4
monoclonal antibody
Product Number: CK64
Aliquot: 1.0 mg
Lot Number: Typical Lot
Storage: -80°C

DESCRIPTION: Monoclonal antibody against Human Bone Morphogenetic Protein 4 (BMP-4)

MOLECULAR WEIGHT: 50 kDa

STORAGE BUFFER: 0.5 x PBS, pH 7.2

STORAGE: Once reconstituted, aliquot and store at -80°C.

IMMUNOGEN: Recombinant Human BMP-4

PURIFICATION: IgG1 antibody was purified by Protein G affinity chromatography.

HOST: Mouse

CLASS: IgG2

APPLICATIONS: ELISA: 1.0 ug/mL
Western: 2.0 ug/mL

SPECIES SPECIFICITY: Human