



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Prothrombin (Mouse) Antiserum

Product Number: CT02

Aliquot: 1.0 mL

Lot Number: Typical Lot

Storage: -70°C

DESCRIPTION:	This polyclonal sheep antibody is specific for mouse prothrombin/thrombin and is provided as antiserum.
MOLECULAR WEIGHT:	N/A
CONCENTRATION:	Whole Serum
SIZE:	1.0 mL
VOLUME:	1.0 mL
STORAGE:	-70°C – Avoid repeated freeze/thaw cycles.
FORM:	Frozen liquid
SOURCE:	Sheep serum
SPECIES SPECIFICITY:	Mouse

For research purposes only.