



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Anti-Thrombopoietin (Mouse) Antiserum

Product Number: CT09

Aliquot: 1.0 mL

Lot Number: Typical Lot

Storage: -70°C

<b>DESCRIPTION:</b>	This polyclonal sheep antibody is specific for mouse thrombopoietin and is provided as antiserum.
<b>MOLECULAR WEIGHT:</b>	N/A
<b>CONCENTRATION:</b>	Whole Serum
<b>SIZE:</b>	1.0 mL
<b>VOLUME:</b>	1.0 mL
<b>STORAGE:</b>	-70°C – <b>Avoid repeated freeze/thaw cycles.</b>
<b>FORM:</b>	Frozen liquid
<b>SOURCE:</b>	Sheep serum
<b>SPECIES SPECIFICITY:</b>	Mouse

For research purposes only.