



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Thrombopoietin (Mouse) FITC
Labeled
Product Number: CT12
Aliquot: 1.0 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	This polyclonal sheep antibody is specific for mouse thrombopoietin and is provided as an FITC labeled IgG fraction.
MOLECULAR WEIGHT:	N/A
CONCENTRATION:	10 mg/mL
SIZE:	1.0 mg
VOLUME:	0.1 mL
STORAGE:	-70°C – Avoid repeated freeze/thaw cycles.
FORM:	Frozen liquid
SOURCE:	Sheep serum
BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
PURIFICATION:	IgG fraction purified by immobilized Protein G.
SPECIES SPECIFICITY:	Mouse

For research purposes only.