

**Product Specification Sheet**

Mouse Prothrombin  
Product Number: CT37  
Aliquot: 0.1 mg  
Lot Number: ct37.z0825  
Storage: -20°C

<b>DESCRIPTION:</b>	Mouse Prothrombin
<b>MOLECULAR WEIGHT:</b>	72 kDa
<b>CONCENTRATION:</b>	8.8 mg/mL
<b>STORAGE BUFFER:</b>	50% Glycerol
<b>STORAGE:</b>	-20°C
<b>PURIFICATION:</b>	Prepared from fresh mouse plasma.
<b>PURITY:</b>	>95% by SDS-PAGE
<b>SPECIFIC ACTIVITY:</b>	3.7 U/mg (1 unitis equivalent to the Factor X activity in 1 mL of normal plasma)