

**Product Specification Sheet**

Description: Sheep Fibrinogen  
Product Number: FB08  
Aliquot: 100 mg  
Lot Number: FB08.750pl  
Storage: -70°C

**DESCRIPTION:** Sheep fibrinogen produced from fresh sheep plasma. Plasminogen was depleted by lysine affinity chromatography.

**MOLECULAR WEIGHT:** 340 kDa

**CONCENTRATION:** Provided lyophilized. Reconstitute With 2.0 mL di water and incubate at 37°C without agitation until completely dissolved. Do not agitate! Mix gently prior to use. Keep aolution at 25-37°C during use. Aliquot and flash freeze unused portion.

**STORAGE BUFFER:** 0.02M Sodium citrate-HCl; pH 7.4

**STORAGE:** -70°C

**SPECIFIC ACTIVITY:** 88% CLOTTABLE

**PURIFICATION:** Purified by multiple chromatographic steps.

**PURITY:** >95% by SDS-PAGE analysis

**APPLICATIONS:**

**SPECIES SPECIFICITY:**