



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Monoclonal Anti-Human Vitronectin
Product Number: GP28
Aliquot: 1.0 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	Monoclonal antibody produced in a vitronectin knockout mouse. Does not cross-react with mouse.
MOLECULAR WEIGHT:	160,000
CONCENTRATION:	10 mg/mL
SIZE:	1.0 mg
VOLUME:	0.1 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Mouse ascites
PURIFICATION:	Immobilized Protein A
APPLICATIONS:	ELISA: Monomeric and multimeric human vitronectin WB: Non-reduced monomeric vitronectin
BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6

For research purposes only.