



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Monoclonal Anti-Human Vitronectin  
Product Number: GP30  
Aliquot: 1.0 mg  
Lot Number: Typical Lot  
Storage: -70°C

<b>DESCRIPTION:</b>	Monoclonal antibody produced in a vitronectin knockout mouse. Does not cross-react with mouse.
<b>MOLECULAR WEIGHT:</b>	160,000
<b>CONCENTRATION:</b>	10 mg/mL
<b>SIZE:</b>	1.0 mg
<b>VOLUME:</b>	0.1 mL
<b>STORAGE:</b>	-70°C
<b>FORM:</b>	Frozen liquid
<b>SOURCE:</b>	Mouse ascites
<b>BUFFER:</b>	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>PURIFICATION:</b>	IgG fraction purified by immobilized Protein A.
<b>APPLICATIONS:</b>	WB: Reduced and non-reduced monomeric vitronectin ELISA: Not suitable.

For research purposes only.