



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA  
USA: 800-692-4633 • Fax: 248-852-4466  
www.oxfordbiomed.com

*Product specifications*

**cDNA: Platelet-Type 12-Lipoxygenase**

**(Murine)**

**Product # MP 12      Typical Lot**

(1.0 µg in 10 µL of TE buffer)

This cDNA can be used for the synthesis,  
by random priming or nick translation,  
of radiolabeled probes that can subsequently  
be used to screen tissues for differential expression  
of 12-LOX isoform

Size of the cDNA = 1.2 Kb

Storage and Stability: At least one year if stored at -70<sup>o</sup> C

Avoid repeated freeze-thaw cycles.

Reference: X-S Chen et.al., JBC 269, 13979-13987, 1994