



OXFORD BIOMEDICAL RESEARCH
P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Product Specification

Single Chain tPA (FITC Labeled)

Product #: PA 52

Lot #: 202

Molecular Weight: 70 KDa

Physical Appearance: yellow frozen solution

Solubility: >2 mg/mL

IgG concentration: 1.44 mg/mL

Considerations:

Storage Conditions: Store frozen at -70°C – stability > 1yr

AVOID REPEATED FREEZE THAW CYCLES

Buffer: 0.4 M HEPES; 0.1 M NaCl; 1 mM EDTA; pH 7.4

Other: >98% pure by SDS PAGE

Fully complexable with human PAI-1