



OXFORD BIOMEDICAL RESEARCH
P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Product Specification

Anti Human Urokinase (uPA)
(rabbit polyclonal)

Product #: PA 65

Lot #: 403

Physical Appearance: colorless frozen solution

IgG concentration: 7.0 mg/mL

Applications:

Western Blotting: 1:100 dilution

Elisa: 1:10,000 dilution

IHC: ~1:10,000 dilution

Considerations:

Storage Conditions: Store frozen at -70°C – stability > 1yr

AVOID REPEATED FREEZE THAW CYCLES