

Product Specification Sheet

Description: Human Urokinase,
HRP labeled
Product Number: PA83
Aliquot: 0.1 mg
Lot Number: pa83.512
Storage: -70°C

DESCRIPTION:	Human urokinase (two chain) that has been conjugated to horse radish peroxidase (HRP).
MOLECULAR WEIGHT:	54,000
CONCENTRATION:	0.2 mg/mL
STORAGE BUFFER:	0.05 M Sodium Acetate; 0.1 M NaCl; 1 mM EDTA; 50% glycerol; pH 5.0
PURITY:	>95% by SDS-PAGE analysis
STORAGE:	-70°C
SOURCE:	Human urine