



## OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA

USA: 800-692-4633 • Fax: 248-852-4466

[www.oxfordbiomed.com](http://www.oxfordbiomed.com)

### Anti-PGH Synthase

**Product # PG 19**      **Typical Lot**

- Immunoglobulin type:** Polyclonal goat IgG.
- Specificity:** Reacts with PGH Synthase from sheep and most other mammalian species including rodents.
- Uses:** Immunoassay, Immunohistochemistry, and Western blots.
- Contents (per vial):** 100  $\mu$ L of goat antiserum containing 0.05% sodium azide.
- Stability:** Stable for several weeks at 4°C, several months at -20°C, and indefinitely at -70°C. Avoid repeated freezing and thawing.
- Protein Concentration:** Approximately 6.0 mg/mL.

**References:**

1. Smith, W.L. and Wilkin, G.P., *Prostaglandins* **13**, 873-892 (1977).
2. Lee, Makau., *Biochemica*, **2**, 19-21, (1996).