



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Human PAI-1 IgG FL Labeled


Product Number: PI 33

Aliquot: 1.0 mg

Lot Number: pi33.1102A

Storage: -70°C

DESCRIPTION:	This rabbit polyclonal antibody is specific for human plasminogen activator inhibitor-1 (PAI-1). The IgG preparation has been fluorescent labeled to increase the utility of the antibody.
MOLECULAR WEIGHT:	50 kDa glycosylated 43 kDa non-glycosylated
CONCENTRATION:	2.5 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.5 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-80°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Recombinant human PAI-1
PURIFICATION:	Purified by Protein A Sepharose.
ISOTYPE:	IgG
APPLICATIONS:	Western: 50-100 µg/mL. ELISA: 1-10 µg/mL.
SPECIES SPECIFICITY:	Human



Western Blot - 10% SDS-PAGE
ASHPAI-GF-BIO-903 at 25 µg/ml Primary
Avidin Alk Phos at 1:100,000 Secondary
1 Human PAI-1 (0.5 µg) Reduced
2 Prestained Standard