



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Anti-Human Complexed Form PAI-1 MoAB IgG

Product Number: PI 63

Aliquot: 1.0 mg

Lot Number: pi63.1104A

Storage: -70°C

<b>DESCRIPTION:</b>	This mouse monoclonal antibody is specific for complexed forms of human plasminogen activator inhibitor-1 (PAI-1). This antibody exhibits neoantigenic properties that allow for binding with any complexed form of PAI-1.
<b>MOLECULAR WEIGHT:</b>	43 kDa non-glycosylated PAI-1 50 kDa glycosylated PAI-1
<b>CONCENTRATION:</b>	2.5 mg/mL as determined spectrophotometrically.
<b>STORAGE BUFFER:</b>	0.1M HEPES; 0.1M NaCl; 1mM EDTA; pH 7.4
<b>STORAGE:</b>	-70°C. <b>AVOID REPEATED FREEZE THAWS.</b>
<b>IMMUNOGEN:</b>	Complexed Recombinant Human PAI-1 and tPA
<b>CLASS:</b>	IgG
<b>CLONE:</b>	MA14D5
<b>PURIFICATION:</b>	Protein A sepharose
<b>APPLICATIONS:</b>	ELISA
<b>SPECIES SPECIFICITY:</b>	Human