



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA

USA 800-692-4633 • Fax 248-852-4466

www.oxfordbiomed.com

Product specifications

Protein Kinase C Gamma (γ) Synthetic Peptide

Product # PK 27 **Typical Lot**

Peptide Type:	Distinct region of protein kinase C C-terminal Sequence: NH ₂ -C-V-H-P-D-A-R-S-P-I-S-P-T-P-V-P-V-M-COOH
Applications:	Will compete with the anti-peptide antibody in Western blotting for confirmation of the PKC specificity of the antibody binding.
Contents per vial:	100 μ g of peptide in 500 μ L of ammonium bicarbonate buffer
Storage:	Stable at 4°C. For long term storage, aliquot and freeze at -20°C. Avoid repeated freeze/thaw cycles.