

Product Specification Sheet

Bovine Plasminogen
Product Number: PL24
Aliquot: 1.0 mg
Lot Number: pl24.411
Storage: -70°C

| | |
|--------------------------|---|
| DESCRIPTION: | Bovine Plasminogen |
| MOLECULAR WEIGHT: | 92,000 |
| CONCENTRATION: | 4.5 mg/mL |
| VOLUME: | 0.222 mL |
| STORAGE BUFFER: | 0.1 M HEPES; 0.1 M NaCl; pH 7.4 |
| STORAGE: | -70°C |
| PURIFICATION: | Prepared from fresh bovine plasma by immobilized lysine chromatography. >95% by SDS-PAGE. |