



OXFORD BIOMEDICAL RESEARCH
P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Product Specification

Rabbit Plasminogen

Product # PL 36

Lot # RbPLG-xxx

Molecular Weight: 92 KDa

Physical Appearance: colorless frozen solution

Solubility: > 2 mg/mL and < 5 mg/mL

Spectrophotometric Data:

Ultraviolet: Absorbance (280nm) = 3.96 $^{0.1\%}$ = 1.69

Active Concentration: 1.89 mg/mL

Activity:

No plasmin activity detected with a chromogenic substrate. > 98% conversion to plasmin observed with activation via uPA.

Considerations:

Storage Conditions: Store frozen at -70°C – stability > 1yr;

Avoid repeated freeze-thaw cycles

Buffer: 0.1 M HEPES; 0.1 M NaCl; 1 mM EDTA; pH 7.4

Other: > 98% pure by SDS-PAGE