

Product Specification Sheet

Description: Human PAI-1,
Glycosylated Latent Form
Product Number: PL77
Aliquot: 0.5 mg
Lot Number: pl77.1104
Storage: -70°C

DESCRIPTION:	Latent Glycosylated wild type human PAI-1 produced in Drosophila S2 cells.
MOLECULAR WEIGHT:	46 kDa
CONCENTRATION:	2.4 mg/ml
STORAGE BUFFER:	0.1M HEPES; 0.1M NaCl; pH 7.4
STORAGE:	Store at -70°C
PURIFICATION:	Separated from active PAI-1 by hydrophobic interaction chromatography.
PURITY:	>95% by SDS-PAGE analysis