



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Product specifications

Rabbit anti-protein serine/threonine phosphatase 1 β

Product # PP 13 Typical Lot

- Description:** Antibody corresponding to a carboxy-terminal peptide of PP1 β catalytic subunit.
- Immunogen:** Purified peptide conjugated to KLH.
Sequence of peptide NH₂-Pro-Arg-Thr-Ala-Asn-Pro-Pro-Lys-Lys-Arg-COOH.
- Specificity:** PP1 β (38 kDa)
- Cross-reactivity:** No cross-reactivity with other recombinant PP1 isoforms.
- Presentation:** 100 μ L whole serum containing 0.02% sodium azide.
- Storage:** Store at -20°C. once thawed store at 4°C.
Avoid repeated freeze-thaw cycles.
Shelf life: 12 months @ -20°C. 6 months @ 4°C.
- Applications:** Immunoblotting - Recommended dilution 1:250 to 1:500
May be used in immunoprecipitation techniques.

For research purposes only