

Product Specification Sheet

Rabbit Anti-Human Prorenin,
Antiserum
Product Number: RN25
Aliquot: 1.0 mL
Lot Number: rn25.310
Storage: -70°C

DESCRIPTION:	This antiserum was produced in rabbit and is against human prorenin. The antiserum recognizes renin and prorenin.
MOLECULAR WEIGHT:	43.7 kDa
STORAGE:	-70°C; minimize freeze/thaw cycles
IMMUNOGEN:	Recombinant full length wt human glycosylated protein
PURITY:	Whole Antiserum
APPLICATIONS:	Western Blot and ELISA
SPECIES SPECIFICITY:	Human