

**Product Specification Sheet**

Description: Rabbit Anti-Human Prorenin Antibody, IgG Fraction, HRP Labeled

Product Number: RN29

Aliquot: 1.0 mg

Lot Number: m29.213

Storage: -70°C

<b>DESCRIPTION:</b>	This rabbit antibody against Human Prorenin has been labeled with HRP and is the IgG fraction. It recognizes Renin and Prorenin.
<b>MOLECULAR WEIGHT:</b>	43.7 kDa
<b>CONCENTRATION:</b>	1.684 mg/mL
<b>STORAGE BUFFER:</b>	PBS, pH 7.4; 1% BSA, 0.01% Thimerosal
<b>STORAGE:</b>	-70°C
<b>IMMUNOGEN:</b>	Recombinant full length wild-type human prorenin.
<b>PURIFICATION:</b>	IgG Fraction.
<b>CLASS:</b>	IgG
<b>APPLICATIONS:</b>	ELISA
<b>SPECIES SPECIFICITY:</b>	Human