



## Product Specification Sheet

Description: Human  
renin/prorenin depleted plasma  
Product Number: RN32  
Aliquot: 10 mL  
Lot Number: m34.1212  
Storage: -70°C

<b>DESCRIPTION:</b>	This human plasma has been doubly depleted of renin and prorenin.
<b>CONCENTRATION:</b>	<0.1 ng/mL Renin, <0.1 ng/mL Prorenin
<b>SOURCE:</b>	Human plasma
<b>STORAGE:</b>	-70°C
<b>PURIFICATION:</b>	Prepared from human plasma using immobilized anti-human renin IgG.