

Product Specification Sheet

Monoclonal antibody against
human rennin and prorenin
Product Number: RN43
Aliquot: 0.1 mg
Lot Number: rn43.509
Storage: -70°C

DESCRIPTION:	This monoclonal antibody was made against human prorenin and binds with both rennin and prorenin.
MOLECULAR WEIGHT:	43.7 kDa
CONCENTRATION:	5.37 mg/mL
STORAGE BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
STORAGE:	-70°C
IMMUNOGEN:	Recombinant full length wt Human glycosylated protein
PURIFICATION:	IgG Fraction
CLASS:	IgG1
CLONE:	7D3-E3
APPLICATIONS:	Western Blotting and ELISA
SPECIES SPECIFICITY:	Has only been tested with Human samples