

**Product Specification Sheet**

Description: Anti Mouse and Rat uPA, IgG Fraction, Biotin Labeled  
Product Number: UP03  
Aliquot: 1.0 mg  
Lot Number: Typical  
Storage: -70°C

<b>DESCRIPTION:</b>	This anti Mouse and Rat polyclonal antibody was produced in Rabbit and is an IgG purification. It has been biotin labeled.
<b>MOLECULAR WEIGHT:</b>	44 kDa
<b>CONCENTRATION:</b>	10 mg/mL
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate, 0.1 M NaCl, 1 mM EDTA, pH 6.6
<b>STORAGE:</b>	-70°C
<b>IMMUNOGEN:</b>	Recombinant full length wild type Mouse glycosylated protein.
<b>PURIFICATION:</b>	Purified by immobilized Protein A.
<b>CLASS:</b>	IgG
<b>APPLICATIONS:</b>	Western Blotting, ELISA
<b>SPECIES SPECIFICITY:</b>	Mouse and Rat