

Product Specification Sheet

Mouse Urokinase LMW
Product Number: UP37
Aliquot: 0.05 mg
Lot Number: up37.1204
Storage: -70°C

DESCRIPTION:	Low molecular weight recombinant active mouse urokinase. Produced in insect cells.
SIZE:	0.05 mg
CONCENTRATION:	1 mg/mL
VOLUME:	0.05 mL
MOLECULAR WEIGHT:	31,000
BUFFER:	0.5M Sodium Acetate; 0.1M NaCl; 1mM EDTA; pH 5.0
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Insect cell culture
PURITY:	>95% by SDS-PAGE analysis

For Research Purposes Only.