

Product Specification Sheet

Description: Mouse PAI-1 tPA
Complex
Product Number: UP57
Aliquot: 0.1 mg
Lot Number: up57.110b
Storage: -70°C

DESCRIPTION: This mouse tPA (single chain) has been reacted with an excess of mouse PAI-1 resulting in a 1:1 molar complex. It has been purified from the remaining PAI-1 by affinity chromatography.

MOLECULAR WEIGHT: 113 kDa

CONCENTRATION: 0.95 mg/mL

STORAGE BUFFER: 0.4M HEPES, 0.1M NaCl, pH 7.4

STORAGE: -70°C

PURIFICATION: Affinity purification.

PURITY: Greater than 95% by SDS-PAGE analysis.