



## Product Specification Sheet

Description: Active Rat  
Urokinase, HMW  
Product Number: VA57  
Aliquot: 0.05 mg  
Lot Number: va57.707  
Storage: -70°C

<b>DESCRIPTION:</b>	This rat urokinase is an active, two-chain, HMW enzyme produced recombinantly in insect cells.
<b>MOLECULAR WEIGHT:</b>	44 kDa
<b>CONCENTRATION:</b>	5.5 mg/ml
<b>STORAGE BUFFER:</b>	0.05M Sodium Acetate; 0.1M NaCl; 1mM EDTA; pH 5.0
<b>STORAGE:</b>	Store at -70°C. Minimize freeze-thaw cycles.
<b>PURITY:</b>	>95% by SDS-PAGE analysis